

# Now Available The All New C·VUE® Advanced™ Custom Toric

**Exceptional Delivery  
Completely Customizable**

**Free Trial Lenses  
Monthly Replacement**

**NEW material with  
Advanced Hydration**



Free Trial Lenses

New Material With Advanced Hydration Properties  
(hioxifilcon D)

Exceptional Deliverability (typically ships in 2 to 3 days)

New Convenient Blister Packaging

Most Customizable Options

Monthly Replacement

+10.00D to -20.00D

Axis in 1° steps

## Another new C·VUE product designed to bring profitability to your practice

AVAILABLE THROUGH LICENSED EYE CARE PROFESSIONALS ONLY

At Unilens, The Independent Eye Care Professional's Specialty Lens Company, our focus is on providing superior contact lens technology while increasing your profitability and enhancing your patient retention. That's why we are now pleased to offer you the C·VUE Advanced Custom Toric in new convenient blister packaging featuring hioxifilcon D, a material which offers advanced hydration properties that allow it to retain 97% of its moisture resulting in outstanding All-Day Comfort. All of these patient benefits, along with exceptional deliverability (lenses typically ship within 2 to 3 days) and our Free Trial/No Returns Required Program helps to contribute to your practice profitability. So call **1-800-446-2020** now to find out how you can advance your practice with the all new C·VUE Advanced Custom Toric!

**Unilens**  
The Independent Eye Care Professional's Lens Company  
www.unilens.com • 800-446-2020

# C·VUE<sup>®</sup> ADVANCED<sup>™</sup> CUSTOM TORIC Soft (hioxifilcon D) Contact Lens Fitting Guidelines

## SUGGESTED PATIENT CRITERIA:

- Normal binocular vision and good ocular health
- Refractive astigmatism between -0.50D and -4.00D

## LENS SELECTION:

### STANDARD PARAMETERS:

Base Curve/Diameter: 8.0/14.0 8.2/14.5 8.5/14.5  
8.8/14.5 9.1/15.0

Powers: +10.00 to -20.00D

Cylinder: up to -4.00D

Axis: any axis in 1° steps

Orientation mark: 6 o'clock position

(Custom parameters available)

### INITIAL FIT DETERMINATION:

#### Base Curve and Diameter

Flattest 'K'	Base Curve	Diameter
47.75D and steeper	8.0	14.0
47.50D to 44.75D	8.2	14.5
44.50D to 42.25D	8.5	14.5
42.00D to 40.00D	8.8	14.5
39.75 and flatter	9.1	15.0

The criteria for the optimum C·VUE Advanced Custom Toric Contact Lens fit is to select the flattest base curve which centers well, exhibits stable rotation and does not move excessively.

#### Power

Order based on the patient's refractive distance Rx in minus cylinder form (vertexed if necessary).

#### Cylinder power

Sum the refractive sphere and cylinder powers and vertex if greater than 4.00D. Calculate the adjusted cylinder power to order by subtracting the vertex corrected distance power.

#### Axis

Order cylinder axis in minus cylinder form.

## LENS FITTING:

- Allow lenses to equilibrate for 10 minutes. Lens comfort should be acceptable after equilibration.
- Lens should center well with 1.0mm to 1.5mm movement with blink in primary gaze.

- Measure acuity in normal room illumination.
- Evaluate the lens rotation. The rotation should be stable from blink to blink.

## SYMPTOM RESOLUTION:

**Excessive movement:** Select steeper base curve.


**Minimal to no movement:** Select flatter base curve.

**Decentration:** Select steeper base curve.

**Inconsistent or excessive rotation:** Select steeper base curve.

**Visual acuity unacceptable:** If rotation is stable, perform a sphero-cylinder over-refraction using the least minus power to obtain acceptable distance vision monocularly. Call Unilens with your sphero-cylinder over-refraction to order the next lens or visit our website at [www.unilens.com](http://www.unilens.com) to use our custom toric lens calculator to determine the correct distance and cylinder powers.

## C·VUE ADVANCED CUSTOM TORIC



*New convenient blister packaging*

**Unilens**

The Independent Eye Care Professional's Lens Company

[www.unilens.com](http://www.unilens.com) • 800-446-2020

TO PLACE AN ORDER OR FOR CONSULTATION  
CALL 1-800-446-2020